



Getting the Best out of Your Automated Liquid Handler

Nathaniel Hentz, PhD

On behalf of Analis Scientific Instruments

LabAutomation: 30 June 2020


Contact: nhentz@artel.co

Today's Agenda

- ❑ Understanding your automated liquid handler (ALH)
- ❑ Linking assay performance to ALH performance
- ❑ ALH measurement and accompanying applications
 - ❑ ALH liquid class optimization
 - ❑ Mixing & plate washing efficiency
 - ❑ Tip evaluation
 - ❑ Training
- ❑ Brief Artel technology description


Why We're Here Today

Global Standard for Liquid Handling Performance



Why We're Here Today

Global Standard for Liquid Handling Performance



Understanding What Automated Liquid Handlers Can Do

Multipurpose automated liquid handling workstations are designed to do any or all assay workflow operations including sampling, heating, cooling, filtering, mixing, washing, diluting and combining of reagents.

- Advantages: minimize contamination, free up lab personnel, reduce assay variability.
- Automated liquid handlers are NOT self-programmable, capable of on-the-fly adjustments, nor are they foolproof.

Automated liquid handling systems must be optimized for each assay

Setting up an Automated Liquid Handler

IQ

Installation Qualification — provides documented evidence and verification that the instrument has been delivered and installed according to manufacturer's specifications

OQ

Operational Qualification — provides documented verification that the instrument subsystems are operating as designed. Verifies that the functionality of an instrument meets the manufacturer's operational specifications.

PQ

Performance Qualification — provides documented verification that the instrument system can perform effectively and reproducibly within performance specifications. Helps ensure confidence in results by verifying that the accuracy and precision of an instrument is maintained.

However...Many Things Affect Liquid Transfer

- Liquid properties (e.g., density, viscosity, surface tension, vapor pressure)
- Liquid transfer device
 - Components (e.g., tips, tubes, reservoirs)
 - Dispense principle (e.g., air/piston displacement, acoustic, pressure, piezo, etc.)
 - Liquid class settings or technique (e.g., air gap, asp/disp speed, immersion depth, etc.)
- Temperature

Performance Capabilities



Tip	Volume (μL)	Inaccuracy (±%)	CV (≤%)
X500	500	1	1
X500	250	3	3
X500	50	5	5
X50	50	1	1
X50	25	3	3
X50	5	5	5

Key to Liquid Class Optimization

- Fluid properties of master mix solutions in various NGS steps are not the same as DNA solutions or diluents!
- Rapid liquid class development for your automated liquid handler
 - Begin by selecting your ALH's default liquid class, such as "PCR", "viscous", "glycerol", etc.
 - Dispense desired volume of test solution
 - Measure dispensed volume
 - Optimize default settings
- Evaluate cold reagent dispense
- Evaluate disposable tips (e.g., lot-to-lot, different vendors)

- pre- and post-air gap volumes
- off-set volume
- aspirate/dispense rate
- aspirate/dispense height
 - tip withdrawal speed
 - wet vs. dry dispense
 - dispense order
 - single dispense or multi-dispense
- tips/cannulas
 - dry or wet tip
 - new or used tip
 - tip touches

Using MVS for Liquid Class Optimization

Default ALH Liquid Class

STATISTICS

Before

Target Volume (µL)	5
Mean volume for all channels (µL)	3.7170
Relative inaccuracy for all channels	-25.66%
Standard deviation for all channels (µL)	0.3907
Coefficient of variation (CV) for all channels	10.51%
Relative inaccuracy pass/fail limit	5%
Coefficient of variation pass/fail limit	5%
Status based on channel results	Failed
Status based on run statistics	Failed

Optimized ALH Liquid Class with MVS®

STATISTICS


After

Target Volume (µL)	5
Mean volume for all channels (µL)	4.9980
Relative inaccuracy for all channels	-0.04%
Standard deviation for all channels (µL)	0.0431
Coefficient of variation (CV) for all channels	0.86%
Relative inaccuracy pass/fail limit	5%
Coefficient of variation pass/fail limit	5%
Status based on channel results	Passed
Status based on run statistics	Passed

Default ALH liquid classes are **not optimal** for all assay solutions, yet it is so **critical**.


MVS Technology Description

- Employs a **dual-dye, dual-wavelength, ratiometric absorbance** measurement for calculating the dispense volume.
- How it works: dyes of known concentration are dispensed into well-characterized microtiter plate. The plate is mixed on a plate shaker to ensure solution homogeneity. Absorbance readings are taken at 520 nm and 730 nm.




Basic Artel MVS Steps

Dispense **red** solution



Add **blue** diluent



Shake on validated MVS plate shaker





Report provides dispense accuracy *and* precision; statistics for entire plate, channel by channel, and dispense order.

ARTEL MVS TEST REPORT		Data Manager 3.3.0.8
Traceable Results*		
Date: 06 Sep 2018		
Time: 4:50:45 PM GMT-5		
Operator: Administrator (admin)		
Liquid Handler Device ID: 8-Channel Device		
Liquid Handler Device Description: Sample 8-Channel Dispensing Tool		
Layout ID: 10uL 8-Channel into 96 Wells		
Layout Description: Dispense 10uL by column into 96-well plate		
Channels: 8		
Plate Description: 96-well MVS Verification Plate		
Dispense Direction: Left to Right		
Device Orientation: Vertical		
Group 1 Statistics		
Target volume (µL)	10	
Target solution	Range B	
Number of data points per channel	12	
Mean volume for all channels (µL)	9.7588	
Relative inaccuracy for all channels	2.41%	
Standard deviation for all channels (µL)	0.3756	
Coefficient of variation (CV) for all channels	3.85%	
Relative inaccuracy pass/fail limit	4%	
Coefficient of variation pass/fail limit	4%	
Status based on channel results	Failed	
Status based on run statistics	Passed	




Read plate

Data Quality is Dependent on Input!



or



Sample/standard preparation

Assay assembly steps

Consumables


















How Do We Employ MVS for an Assay?

Typical ELISA protocol:

1. Coat plate with 200 μ L of Capture antibody, incubate overnight at 4°C
2. Wash 4x with 200 μ L of Wash buffer
3. Block plate with 300 μ L of Block solution for 2 h at 37°C
4. Wash 4x with 300 μ L of Wash buffer
5. Add 100 μ L of sample, prediluted 1:100 in Block solution
 - a. Seal plate
 - b. Incubate 1 h at 37°C
6. Wash 4x with 300 μ L of Wash buffer
7. Dilute Primary antibody 1:100 with Block solution
 - a. Add 100 μ L Primary antibody per well
 - b. Incubate 1 h at 37°C
8. Wash 4x with 300 μ L of Wash buffer
9. Dilute Detection antibody 1:1000 with Block solution
 - a. Add 100 μ L of Secondary antibody per well
 - b. Incubate 30 min at 37°C
10. Wash 4x with 300 μ L of Wash buffer
11. Add 100 μ L of TMB Substrate per well
 - a. Incubate 10 min at 37°C
12. Add 100 μ L of Stop solution per well
13. Read on plate reader at absorbance 450 nm

How Do We Employ MVS for an Assay?

Typical ELISA protocol:

1. Coat plate with 200 μ L of Capture antibody, incubate overnight at 4°C
2. Wash 4x with 200 μ L of Wash buffer 
3. Block plate with 300 μ L of Block solution for 2 h at 37°C 
4. Wash 4x with 300 μ L of Wash buffer 
5. Add 100 μ L of sample, prediluted 1:100 in Block solution   
 - a. Seal plate
 - b. Incubate 1 h at 37°C
6. Wash 4x with 300 μ L of Wash buffer 
7. Dilute Primary antibody 1:100 with Block solution
 - a. Add 100 μ L Primary antibody per well 
 - b. Incubate 1 h at 37°C
8. Wash 4x with 300 μ L of Wash buffer 
9. Dilute Detection antibody 1:1000 with Block solution 
 - a. Add 100 μ L of Secondary antibody per well 
 - b. Incubate 30 min at 37°C
10. Wash 4x with 300 μ L of Wash buffer 
11. Add 100 μ L of TMB Substrate per well  
 - a. Incubate 10 min at 37°C
12. Add 100 μ L of Stop solution per well   
13. Read on plate reader at absorbance 450 nm

Wash verification



Volume verification

Serial dilution accuracy

Mixing verification

Liquid class optimization

Optimizing the Assay Workflow



These are statistically more distinguishable from each other, allowing better decision making

CMPD1, 2, 3, and 4 are structurally related compounds with similar potencies (n=4)

Thank You!

Let's take some questions now

Contact: nhentz@artel.co